

**MAGNETIC STIRRER MS-1 / MS-5 model – ZELTEC -
Made in Argentina**

The line of economic Zeltec magnetic stirrer presents these models capables of shaking from a recipient of 10 ml to 3 liters of a glass.

Engines with protection against spills and vapors do not require maintenance, thus ensuring, together with its modern electronics a long life.

The speed can be adjusted between 60 and 1200 rpm **(1)** in electronic form and can function continuously for hours **(2)**

References

(1) Verified stirring at room temperature water in a recipient of 1 liter.

(2) In both models there was good for 48 hours continuous operation at 600 rpm, room temperature (ms-1) and 80 ° C (ms-5), stirring liquid refrigerant in a glass of 1 liter.



Fitted with thermal fuse self

MS 1 model without thermostatisation

MS 5 model with thermostatisation



FEATURES TECHNIQUES

| | MS 1 | MS 5 | Comments |
|--|----------------------------|-------------------------------|--------------------|
| Function | Magnetic stirrer | Magnetic stirrer with heating | - |
| Alimentation | 110 v 60 Hz 220 v 50 Hz | 110 v 60 Hz 220 v 50 Hz | For export |
| Regulation of speed in rpm | 60 – 1200 | 60 – 1200 | A request may vary |
| Thermostation of temperature | - | Ambient to 300 ° C | - |
| Control Temperature | - | Manual thermostats | - |
| Stainless Steel surface agitation | 145 x 180 mm | 145 x 170 mm | - |
| Power | 33 W | 550 W | - |
| Measures long wide high | 140 mm 180 mm 93 mm | 154 mm 183 mm 118 mm | - |
| Measures Packing long wide high | 270 mm 200 mm 160 mm | 270 mm 200 mm 160 mm | - |
| Weight | 1.8 kg | 2.3 kg | - |

