

**TZC 100 / ZC 200 model – ZELTEC -
Made in Argentina**

This equipment allows thermostatic Zeltec family quickly and accurately a large quantity of tubes. Digital Timer has a program that is starting to work. At met cycle audible alarm sounds a warning and automatically shuts off the power, regardless of the presence of an operator. It comes with three preset temperatures over a manual control.



These models have PT100 sensor, which allows a sophisticated electronics with a power accuracy and stability of ± 0.2 °C in temperature and pre ± 0.5 degrees C in temperature set manually.

The digital timer can be programmed from 1 minute to 19 hours, doing the countdown. 5 minutes later than the time it sounds a warning alarm, and at the end of the cycle longer sounds an alarm, automatically cutting the power.

In the middle of the team, however, there is a recess to accommodate a vertical thermometer.

A thermal fuse protection provides greater security.

The ZC 200 model is compatible with Hach COD tubes, and NANOCOLOR MERCK.





TECHNICAL CHARACTERISTICS

	ZC 100	ZC 200
Temperature Range	Ambient 100 °C	Ambient 200 °C
Preset temperatures	3 temperature choice	100 – 148 – 160 °C
Alimentation	220 V / 50 Hz	220 V / 50 Hz
Power	600 W	700 W
Stream	2,75 Amp	3,18 Amp
measures	Length 230 mm Width 260 mm High 97 mm	Length 230 Width 260 Hight 97
Packing	Length 320 Width 320 Hight 220	Lenght 320 Width 320 Hight 220

